




**STANDS FIRM
IN EVERY SITUATION**



**STABLE RESULTS—
VISALYS® CORE**

**DUAL-CURING COMPOSITE FOR CORE BUILD-UP
AND ROOT POST CEMENTATION**



KETTENBACHDENTAL
Simply intelligent

VISALYS® CORE – EASY TREATMENT

VISALYS® CORE – STABLE RESULTS



Dual-curing core build-up & root post cementation material

Visalys® Core – especially for difficult and freehand core build-ups.



EXCELLENT STABILITY

- with very good flow properties at the same time



UNIQUE ACTIVE-CONNECT-TECHNOLOGY

- for optimized adhesive bonding
- for use with common adhesives



DENTIN-LIKE GRINDABILITY




FREE OF BISPHENOL A




Distinguished by the experts “The Dental Advisor”

Secure bonding with all relevant adhesives



Visalys® Core achieves the clinically relevant adhesive bond of approx. 10 MPa.

Removal depth when grinding various core build-up materials



Visalys® Core can be optimally processed due to its dentin-like grindability.

VISALYS® CORE – THE RANGE

Visalys® Core

Intro pack	White Dentin	1 x 9 g (5 mL) automix syringe, 10 grey mixing tips, 5 intraoral tips, 5 endo tips	REF 13866 REF 13865
Normal pack	White Dentin	2 x 9 g (5 mL) automix syringe, 20 grey mixing tips, 10 intraoral tips, 10 endo tips	REF 13860 REF 13861
Normal pack	White Dentin	1 x 45 g (25 mL) cartridge, 20 yellow mixing tips, 20 intraoral tips	REF 13870 REF 13871



Kettenbach GmbH & Co. KG
 Im Heerfeld 7
 35713 Eschenburg · Germany
 Telefon +49 (0)2774 7050
 info@kettenbach.com
 www.kettenbach-dental.com



027740_3425


KETTENBACHDENTAL
 Simply intelligent