




EFFICIENT IMPRESSIONS – **SILGINAT®**

A-SILICONE FOR SITUATION IMPRESSIONS



KETTENBACHDENTAL

Simply intelligent

SILGINAT® – SIMPLE, FAST, EFFICIENT!

Anatomical impression

Silginat® is an alginate replacement material based on A-silicone especially for **alginate applications**.



INDICATIONS

- Opposing dentitions
- Fabricating temporary crowns and bridges
- Anatomical models
- Fabricating simple removable prosthetic restorations
- Producing models for the construction of trays for whitening solutions, bites and surgical guides
- Orthodontic work
- Case study models
- Fabricating removable retainers and splints



Technical data



Working time
at 23 °C /
74 °F ≤

Intraoral
setting time at
35°C / 95 °F ≥

Total
setting
time* ≥

Elastic
recovery test
(approx.) %

Linear
dimensional
change
(approx.) %

Hardness (ap-
prox.) Shore



Silginat® /
Silginat® Strawberry
1:1 cartridge

1:15 min.

1:15 min.

2:30 min.

99.5

– 0.20

A 42



Silginat® /
Silginat® Strawberry
5:1 cartridge

1:30 min.

1:30 min.

3:00 min.

99.5

– 0.20

A 45

*Total setting time (removal time from mouth) from the start of mixing

AN A-SILICONE AS YOU KNOW IT



- Storage stability and dimensional stability, high precision

Silginat vs. alginate
approx. 24 h after impression-taking



EASY TO HANDLE



- Clean and fast work thanks to the 5:1 jumbo or 1:1 50 mL cartridge with reproducible results



COST-EFFECTIVE IMPRESSION TAKING

- Pourable multiple times, only one anatomical impression per patient required



FOR EVEN GREATER EFFICIENCY

Saves time and many work steps during your work!

SILGINAT® – THE RANGE

Silginat®		
Refill pack	2 x 380 mL cartridges	REF 14713
Normal pack	6 x 50 mL cartridges, 6 mixing tips	REF 13846
Bulk pack	24 x 50 mL cartridges	REF 13847

Silginat® Strawberry		
Refill pack	2 x 380 mL cartridges	REF 14715
Normal pack	6 x 50 mL cartridges, 6 mixing tips	REF 13826
Bulk pack	24 x 50 mL cartridges	REF 13827

Accessories		
Mixing tips white, Ø 6,5 mm	50 tips	REF 17234
Mixing tips white, Ø 6,5 mm	100 tips	REF 17235
Dynamic mixers white	45 tips	REF 17900



Kettenbach GmbH & Co. KG
Im Heerfeld 7
35713 Eschenburg · Germany
Phone +49 (0)2774 7050
info@kettenbach.com
www.kettenbach-dental.com



1410873_3625

KETTENBACHDENTAL
Simply intelligent