



SOOO COMFORTABLE IS **IDENTIUM**[®]

PRECISION IMPRESSION



KETTENBACHDENTAL

Simply intelligent

IDENTIUM[®] – THE UNIQUE MATERIAL COMBINATION!

A new generation of precision impression material

Combines the advantages of two proven impression materials (A-silicone and polyether) – Vinylsiloxanether[®].




ONE-STEP IMPRESSION

- Double-mix impressions
- Monophase impressions



INDICATIONS

- Crowns and bridges
- Inlays and onlays
- Functional impressions
- Implants
- Combined restorations (fixation impressions)
- Pick-up impressions



Double-mix impression




Implant impression



IDENTIUM[®] – COMFORTABLE IN ALL DIMENSIONS

Removal from the mouth - as light as a feather

Identium[®] vs. Impregum[™]: **2.6 times** lower withdrawal forces



Ease of removal from the patient's mouth avoids distortions and ensures a stress-free impression taking process.


For your precision.



Advantages your patient will appreciate too

Fast final setting time

with sufficient processing time and short intraoral setting time



Sufficient intraoral working time makes it ideal for large cases.
For your precision.

Outstanding flowability



With the outstanding flowability of **Identium® Light**, you can reproduce even the narrowest gaps in great detail. **For your precision.**



- **Odorless and tasteless material** reduces gag reflex
- **Quick setting and easy removal from the mouth** for a relaxed and pleasant experience
- **Initially hydrophilic** – good performance even with residual moisture

COMFORTABLE AND PRECISE IMPRESSIONS



EASY REMOVAL FROM THE MOUTH

SHORT INTRAORAL SETTING TIME

INITIALLY HYDROPHILIC

**Immediately start
applying the material:**
no discarding, no waste



Available in viscosities
Heavy, Medium and Light



Selectable setting times
regular (4:30 min.)
fast (3:30 min.)



Compatible with most common
mixing devices



IDENTIUM® – THE RANGE

		Identium®			
HEAVY		Identium® Heavy	Identium® Heavy Fast		
	Intro pack	REF 1472411	–		
		380 mL cartridge, 50 mL cartridge Identium® Light, 10 mL adhesive, 10 dynamic mixers, 8 red mixing tips, 8 transparent intraoral tips	–		
	Refill pack	REF 14725	REF 14727		
2 × 380 mL cartridges		2 × 380 mL cartridges			
MEDIUM		Identium® Medium	Identium® Medium Fast		
	Intro pack	REF 14716	–		
		380 mL cartridge, 10 mL adhesive, 10 dynamic mixers, 1 application syringe	–		
	Refill pack	REF 14717	REF 14719		
2 × 380 mL cartridges		2 × 380 mL cartridges			
LIGHT		Identium® Light	Identium® Light Fast		
	Normal pack	REF 1370111	REF 1371111		
2 × 50 mL cartridges, 8 red mixing tips		2 × 50 mL cartridges, 8 red mixing tips			
ACCESSORIES	Identium® Adhesive	Dynamic Mixers	Mixing tips	Sympress Dispenser	Applyfix 5
	REF 14204	REF 17900	REF 17249 / REF 17250	REF 35910	REF 17204
	10 ml bottle	45 pieces blue	50 mixing tips red / 100 mixing tips red	Electrical dosing and mixing device	2 application syringes, 12 syringe tips



NOTES

Kettenbach SNC
7, Place de la Gare
57200 Sarreguemines · France
info@kettenbach.com
www.kettenbach-dental.com

Not for use, sale and
distribution in Germany.



KETTENBACHDENTAL
Simply intelligent



1410870_1024