

THE LUTING COMPOSITE FOR
ALL INDIRECT
RESTORATIONS



SOOO SMART
IS **VISALYS**® CEMCORE

CEMENTATION & CORE BUILD-UP



KETTENBACHDENTAL
Simply intelligent

VISALYS® CEMCORE – ENTIRELY CONVINCING

FOR ALL CASES

- **Suitable for permanent cementation** of crowns, bridges, inlays, onlays, veneers and adhesive bridges, root posts and indirect core build-ups
- **Suitable for all dental restoration materials** ((noble)metal, oxide ceramic, (zirconium oxide and aluminum oxide ceramics), silicate ceramic (feldspar and glass ceramic), hybrid ceramic, composite)

Conclusion by Prof. Dr. Claus-Peter Ernst

[...] (One) impressive potential by **Visalys® CemCore** for a huge range of indications [...]

Source: Ernst, Prof. Dr. Claus-Peter et al.: A universal adhesive cementation and core build-up material, in: Spitta: ZMK, Vol. 36, Issue 5 May 2020 [294–308]



HIGH COLOR STABILITY

Unlike other luting composites, **Visalys® CemCore** is hydrophobic. This minimizes the absorption of colors and discoloration of the cement joint.



	RelyX™ Universal	RelyX™ Unicem	RelyX™ Ultimate	Panavia™ F2.0	Panavia™ V5	Multilink® Automix	Variolink® Esthetic DC	Visalys® CemCore
Water								
Coffee								
Red wine								

Light-cured samples, stored 48 h in the corresponding solution at 37 °C.

SECURE ADHESION & LONG LIFE





The **high quality of the marginal gap** of Visalys® CemCore represents the high quality of the bond between the dental substance and the composite.

Source: Prof. Dr. Roland Frankenberger et al. (2019), Marginal Quality of Ceramic Inlays luted with Visalys® CemCore, University Marburg

No debonding and no fracture, even after **10 years** of simulated oral use.

Source: Prof. Dr. Dipl.-Ing. Martin Rosentritt (2020), Fracture force after long-term storage and TCML, University Regensburg



SMART HANDLING

Always the same procedure, regardless of the type of indirect restoration.



Easy to remove any excess

1) Outstanding adhesion results with patented Active-Connect-Technology.

More Information: www.kettenbach-dental.com/act

2) Copyright: Dr. med. dent. Frank-Michael Maier



Now view the video
"Fixation of 16 veneers"

VISALYS® CEMCORE – THE RANGE



Visalys® CemCore

Starter pack

4.5 g (2.5 mL) automix syringe,
2 mL Tooth Primer, 2 mL Restorative Primer,
blue mixing tips, blunt/tapered, 5 units each,
3 intraoral tips, 2 endo tips

Universal (A2/A3)

REF 13570

Normal pack

9 g (5 mL) automix syringe,
blue mixing tips, blunt/tapered, 10 units each,
6 intraoral tips, 4 endo tips

Universal (A2/A3)

Opaque

REF 13572

REF 13573

Translucent

REF 13574

Bleach

REF 13575

Dark (A4)

REF 13576

Visalys® Tooth Primer

4 ml bottle

REF 13580

Visalys® Restorative Primer

4 ml bottle

REF 13581

Visalys® CemCore Try In Paste

1.4 mL syringe,
5 application tips

Universal (A2/A3)

Opaque

REF 13592

REF 13593

Translucent


REF 13590

Bleach

REF 13591

Dark (A4)

REF 13594



Made
in
Germany

NOTES



KETTENBACHDENTAL

Simply intelligent

Kettenbach SNC
7, Place de la Gare
57200 Sarreguemines · France
info@kettenbach.com
www.kettenbach-dental.com

1410524_1124

