



SOOO STABLE
IS **VISALYS[®]** CORE

CORE BUILD-UP & POST CEMENTATION



KETTENBACHDENTAL

Simply intelligent

VISALYS® CORE – EASY TREATMENT

Dual-curing core build-up & root post cementation material

Visalys® Core – especially for difficult and freehand core build-ups.



EXCELLENT STABILITY

- with very good flow properties at the same time



UNIQUE ACTIVE-CONNECT-TECHNOLOGY

- for optimized adhesive bonding
- for use with common adhesives



DENTIN-LIKE GRINDABILITY




FREE OF BISPHENOL A



Distinguished by the
experts “The Dental
Advisor”


VISALYS® CORE – STABLE RESULTS

Secure bonding with all relevant adhesives



Visalys® Core achieves the clinically relevant adhesive bond of approx. 10 MPa.

Removal depth when grinding various core build-up materials



Visalys® Core can be precisely processed similar to dentin.

VISALYS® CORE – THE RANGE



Visalys® Core

Intro pack	White Dentin	1 x 9 g (5 mL) automix syringe, 10 brown mixing tips, 5 intraoral tips, 5 endo tips
Normal pack	White Dentin	2 x 9 g (5 mL) automix syringe, 20 brown mixing tips, 10 intraoral tips, 10 endo tips
Normal pack	White Dentin	1 x 45 g (25 mL) cartridge, 20 yellow mixing tips, 20 intraoral tips



NOTES

Kettenbach SNC
7, Place de la Gare
57200 Sarreguemines · France
info@kettenbach.com
www.kettenbach-dental.com



022740_1024


KETTENBACHDENTAL
Simply intelligent