



# HIGHLY IMPACTFUL – **VISALYS® TEMP**

**TEMPORARY CROWN AND BRIDGE MATERIAL**



**KETTENBACHDENTAL**

Simply intelligent

# VISALYS® TEMP – ABSOLUTELY BRILLIANT

## Temporary crown and bridge material

Visalys® Temp – for the fabrication of temporary crowns, partial crowns, bridges, inlays, onlays and veneers.

HIGHLY DURABLE

ALSO SUITABLE AS LONG-TERM  
TEMPORARIES

HIGHLY POLISHABLE

EASY & PRECISE PROCESSING


EASY TO REPAIR



Distinguished by the experts  
"The Dental Advisor"

# VISALYS® TEMP – HIGHLY IMPACTFUL


## High breaking strength



Measurements from the University Hospital Regensburg, based on EN ISO 20795-1:2013 (D).

High breaking strength reduces the risk of a temporary denture breaking in the mouth under high or sudden chewing loads

## Working and setting times



**Note:**

Working times apply for products, which have been stored and processed at a temperature of 23°C/74°F and a relative humidity of 50%. Higher temperatures shorten the working time and lower temperatures extend these times.

# VISALYS® TEMP – THE RANGE

## Visalys® Temp

<b>Normal pack</b>	A1		REF 13780
	A2		REF 13781
	A3	1 x 71 g (50 mL) cartridge, 15 blue-orange mixing tips	REF 13782
	A3.5		REF 13790
	B1		REF 13784
	BL		REF 13788
<b>Bonus pack</b>	A2	5 x 71 g x (50 mL) cartridges, 15 blue-orange mixing tips	REF 13794
	A3		REF 13795



Made  
in  
Germany

142741\_4825



Kettenbach GmbH & Co. KG  
Im Heerfeld 7  
35713 Eschenburg · Germany  
Phone +49 (0)2774 7050  
info@kettenbach.com  
www.kettenbach-dental.com

  
**KETTENBACHDENTAL**  
Simply intelligent